

EcoCRM™

A Recombinant CRM₁₉₇ Carrier Protein

A. Lees¹, R. Simon², S. Baliban², I. Krauss³, D. Nguyen³, M. Pravetoni^{4,5}, N. Oganessian¹


^{1,2}Fina Biosolutions, Rockville, MD; ²U Maryland School of Med, Baltimore; ³Brandeis U, Waltham, MA; ⁴Minneapolis Medical Research Foundation (MMRF); ⁵Center for Immunology, U Minnesota

EcoCRM™ Affordable CRM₁₉₇ Carrier Protein

CRM₁₉₇, a genetically detoxified diphtheria toxin, is widely used as a carrier protein in conjugate vaccines. It was originally expressed as a secreted protein in *Corynebacterium diphtheriae*. Until recently, both “native” and recombinant CRM₁₉₇ has been difficult to obtain and/or expensive.

Fina BioSolutions has developed a new CRM₁₉₇, EcoCRM™, produced in *E. coli*. EcoCRM™ is expressed as soluble protein in the cytoplasm of a widely-used *E. coli* expression strain. High expression of EcoCRM™ in a BL21 strain and a simple purification method allows low cost production and the promise of significantly reducing the cost of this component of conjugate vaccines.

The data presented show that EcoCRM™ is an excellent carrier protein for a variety of antigens.




SDS PAGE EcoCRM™
>99% purity. No nicking evident.

SEC HPLC of EcoCRM™
Dimer is <0.6%. Purity >99%.

Therapeutic Vaccines for Treatment of Opioid Abuse and Overdose


The Pravetoni lab has developed a series of conjugate vaccines that elicit opioid-specific IgG antibodies that reduce opioid distribution to the brain. Immunization selectively blocks opioid-induced behaviors, such as opioid self-administration, and prevents opioid-induced respiratory depression and bradycardia in mice and rats. Current efforts focus on readying for FDA approval and Phase I clinical trial a vaccine consisting of oxycodone (OXY) conjugated to a GMP monomer KLH carrier protein. Due to its mullosk origin and poor characterization, native KLH is not an ideal carrier. This study compared oxycodone conjugates using EcoCRM™ and KLH.



The haptenization ratio of oxycodone-EcoCRM™ can be characterized by MALDI-TOF.

EcoCRM™ was analyzed before and after conjugation to the OXY hapten. A ratio of 24 haptens per CRM was calculated from the MW difference. We were unable to characterize the conjugates containing native KLH.

IN VIVO VACCINE EFFICACY



The oxycodone-EcoCRM™ conjugate showed superior efficacy than the previously established oxycodone-KLH. In contrast to KLH conjugates, EcoCRM™ conjugates are easier to characterize.

The Pravetoni group is currently developing vaccines and monoclonal antibodies against heroin, oxycodone, and fentanyl abuse and overdose (Pubmed, “Pravetoni M”). For collaborations or partnership, please contact Dr. Marco Pravetoni at prave001@umn.edu.

Induction of oxycodone-specific IgG antibody prevents oxycodone distribution to the brain in mice.

BALB/c mice (n ≥ 4/group) were immunized s.c. with 100 µg unconjugated KLH (control), OXY-KLH, or OXY-EcoCRM™, formulated with 1 mg of alum adjuvant on days 0, 14 and 28. On day 35, mice were challenged with 2.25 mg/kg oxycodone, and then brain and serum collected for analysis of vaccine efficacy.

A) Oxycodone-specific serum IgG titers analyzed by ELISA.
B) Oxycodone concentration in the brain 30-min after oxycodone challenge, analyzed by GC/MS. Above bars, percentages (%) indicate decrease in brain oxycodone compared to KLH.


Induction of oxycodone-specific B cells

BALB/c mice (n ≥ 4/group) were immunized i.m. with 100 µg KLH, OXY-KLH, or OXY-EcoCRM™, formulated with 200 µg of alum. At 7 days post-immunization, OXY-specific B cell populations from spleen and lymph nodes were isolated by antigen-based magnetic enrichment paired with multi-parameter flow cytometry.

C) Oxycodone-specific IgM+ B cells.
D) Oxycodone-specific germinal center (GC) B cells.

S. Typhimurium Core-OPS-EcoCRM™ Conjugate

Development of a glycoconjugate vaccine to prevent invasive Salmonella Typhimurium infections in sub-Saharan Africa. Baliban et al. PLoS Negl Trop Dis 11(4):e0005493, 2017



Serum IgG titers against STm-COPS from mice (n = 20/group) immunized with 2.5 µg STm-COPS-KDO:CRM₁₉₇ or PBS. Solid bars indicate the GMT.


Kaplan-Meier survival curves of mice immunized with STm-COPS-KDO-CRM₁₉₇ (circle) or PBS (diamond) after challenge (n = 20/group) with 1 or 5x10⁶ CFU of STm D65.

S. Typhimurium COPS-EcoCRM™ conjugate provided 100% protection in this animal model


HIV Glycopeptide-EcoCRM™ Conjugate

Synthesis of multivalent glycopeptide conjugates that mimic an HIV epitope. Bailey et al. Tetrahedron. 72:6091, 2016.

Thio-ether conjugation is frequently used for linking peptides, oligosaccharides and other antigens to a carrier protein such as CRM₁₉₇.




EcoCRM™ was labeled with an excess of the NHS maleimide reagent sulfo-EMCS. The MW of EcoCRM™ increased from 58,413 to 62,525, (delta=4100). Each linker adds a mass of 208, indicating about 20 of CRM's 39 lysines have been labeled.



Thiol-glycopeptide was conjugated to the EMCS-labeled CRM₁₉₇. Mass spec of the conjugates indicated ~7 glycopeptides/CRM.

Rabbit responses to different CRM conjugates



Glycopeptide conjugates of EcoCRM™ and CRM from *P. fluorescens* (“Pfenex”) were prepared. Rabbits were immunized twice. After the final bleed, IgG anti-glycopeptide titers were measured using glycopeptide-BSA as the ELISA antigen.

Glycopeptide conjugates with EcoCRM™ & Pfenex CRM induced comparable anti-glycopeptide titers in rabbits.



Fina Biosolutions
www.FinaBio.com
Info@FinaBio.com

EcoCRM™ is an effective carrier protein for conjugate vaccines and is an alternative to CRM₁₉₇ from *Corynebacterium* and *Pseudomonas*. A simple expression system, high yields and the efficient purification of EcoCRM™ can help reduce the cost of conjugate vaccines.