01

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


5e5 of the Anti-CD19 CAR-293 cells were stained with 100  $\mu$ L of 1:50 dilution (2  $\mu$ L stock solution in 100  $\mu$ L FACS buffer) of PE-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. No. **FM3-HPY53**) and negative control antibody (Cat. No. **DNP-PM1**) respectively, PE signal was used to evaluate the binding activity (QC tested).

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
PE-Labeled Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Cat. No. **FM3-HPY53**) at 25°C for 48 hours and freeze-throw cycles are stable without performance reduction.

5e5 of PBMCs were stained with a series of concentrations of PE-Labeled Monoclonal Mouse IgG1 Isotype Control (Cat. No. **DNP-PM1**), and PE signal was used to evaluate the binding activity (QC tested).

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Binding activity of two different lots of PE-Labeled Monoclonal Anti FMC63 scFv Antibody, Mouse IgG1(Cat. No. **FM3-HPY53**) against Anti-CD19 CAR-293 cells was evaluated by flow cytometry. The result shows very high batch-to-batch consistency.

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


Detection of FMC63 ADA titer by bridging-ELISA Assay. Immobilized Mouse FMC63 scFv can bind FMC63 ADA in 1:10 human serum, and then add Biotin-Mouse FMC63 scFv. Detection was performed using HRP-conjugated streptavidin with sensitivity of 78 ng/mL (QC tested).

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| アッセイの種類            | カタログ番号      | 種      | 説明                                                                                                         |
|--------------------|-------------|--------|------------------------------------------------------------------------------------------------------------|
| ADA Assay          | FM3-S93     | Rabbit | Rabbit Anti-Mouse FMC63 scFv Polyclonal Antibody (recommended for ADA assay) (MALS verified)               |
|                    | CD9-M52Hb   | Mouse  | FMC63 scFv, His Tag (SPR verified)                                                                         |
|                    | CD9-BV4324b | Mouse  | Biotinylated Anti-CD19 scFv, Fc, Avitag™ (FMC63) (MALS verified)                                           |
| Neutralizing Assay | FM3-Y45A1   | Mouse  | Monoclonal Anti-FMC63 scFv Antibody, Mouse IgG1 (Y45) (Carrier-free) (recommended for ADA assay) DMF Filed |

**Acro** BIOSYSTEMS  
+ 1 800-810-0816 (US / Canada)  
+ 86 400-682-2521 (Asia & Pacific)  
[techsupport@acrobiosystems.com](mailto:techsupport@acrobiosystems.com)  
[www.acrobiosystems.com](http://www.acrobiosystems.com)  
1 Innovation Way,  
Newark, DE 19711

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